

# NorthEast Monitoring, Inc.

*advancing Holter technology™*

## Holter LX<sup>®</sup> Analysis Software - Features Matrix

Software Feature / Arrhythmia	Software Level				
	PRO	ENHANCED+	ENHANCED	BASIC	REMOTE
Remote Send,	✓***	✓***	✓***	✓***	✓***
Remote Send, Other	✓	✓	✓	✓	✓
1 Strip Per Hour	✓	✓	✓	✓	
Bi-Directional Scanning	✓	✓	✓	✓	
Diary Entry	✓	✓	✓	✓	
Event Marker	✓	✓	✓	✓	
Label Strips	✓	✓	✓	✓limited	
Maximum Heart Rate	✓	✓	✓	✓	
Minimum Heart Rate	✓	✓	✓	✓	
Page Mode	✓	✓	✓	✓	
SpO2	✓**	✓**	✓**	✓**	
<b>Configuration</b>					
Calibration Strips	✓	✓	✓		
Configuration	✓	✓	✓		
Critical Events	✓	✓	✓		
Failure To Capture	✓	✓	✓		
Failure To Sense	✓	✓	✓		
Inhibition	✓	✓	✓		
Insert Beat	✓	✓	✓		
Interval Editor	✓	✓	✓		
Pause	✓	✓	✓		
Saved Strips	✓	✓	✓		
ST Baseline	✓	✓	✓		
ST Depression	✓	✓	✓		
ST Deviation (Maximum)	✓	✓	✓		
ST Elevation	✓	✓	✓		
SVPB, SVPB Pair, SVT	✓	✓	✓		
Time Domain Hrv	✓	✓	✓		
True 12-Lead	✓*	✓*	✓*		
VPB, VPB Pair, VTAC	✓	✓	✓		
<b>Analysis</b>					
AFib On/Off	✓	✓			
Remote Receive,	✓***	✓***			
Remote Receive, Other	✓	✓			
Review Bins	✓	✓			
Templates	✓	✓			
Beats	✓	✓			
<b>Advanced Analysis</b>					
Bigeminy	✓				
Bradycardia	✓				
Early VPB	✓				
Fastest SVT	✓				
Fastest VTAC	✓				
Irregular R-R	✓				
Longest R-R	✓				
Longest SVT	✓				
Longest VTAC	✓				
Quadragesiminy	✓				
R On T VPB	✓				
Spectral Analysis	✓				
Spectral Plots & Tables	✓				
ST End	✓				
ST Maximum HR	✓				
ST Onset	✓				
Tachycardia	✓				
Trigeminy	✓				

\* DR180+ Recorder required \*\*DR180+ and purchase of SpO2 cable required \*\*\* MagicVORTEX Account required



## REPORT Features Matrix


Report Feature	Software Level			
	PRO	ENHANCED+	ENHANCED	BASIC
Back-up and restore (CD-RW required)	✓	✓	✓	✓
Comments page	✓	✓	✓	✓
Full disclosure print - 1 hour to 48 hours	✓	✓	✓	✓
Full disclosure print - single channel	✓	✓	✓	✓
Full disclosure print - two channel	✓	✓	✓	✓
List of diary events	✓	✓	✓	✓
Patient information summary: 4 choices	✓	✓	✓	1 only
Print full-sized strips	✓	✓	✓	✓
Print review screen	✓	✓	✓	✓
Status window	✓	✓	✓	✓
12-lead table (25 pages)	✓*	✓*	✓*	✓*
12- lead strips - 25 mm/sec	✓*	✓*	✓*	
12-lead strips - 50 mm/sec	✓*	✓*	✓*	
General profile/trends	✓	✓	✓	
HRV time summary	✓	✓	✓	
ST episodes	✓	✓	✓	
Time domain	✓	✓	✓	
Normal bins	✓	✓		
Paced bins	✓	✓		
Paced data information / heart rate trend	✓	✓		
Ventricular bins	✓	✓		
12-lead trends	✓*			
3-D plot "waterfall graph" (spectral)	✓			
Bradycardia/Tachycardia Table	✓			
Circadian plots (spectral)	✓			
Critical events	✓			
Frequency domain table (spectral)	✓			
Hourly rhythm page	✓			
List of saved events	✓			
Normalized frequent domain (spectral)	✓			
Paced interval histograms	✓			
Paced summary	✓			
Print expanded HR/ST trends	✓			
Scanning criteria/settings	✓			
Spectral HRV (frequency domain)	✓			
ST analysis triaxial graph	✓*			
Supraventricular summary	✓			
Ventricular summary	✓			

\*12-lead recorder option required

### Industry-leading 3-year warranty

NB: In the U.S.A., Federal Law restricts devices to sale by or on the order of a physician. For in vitro diagnostic use. Not intended to replace real-time telemetry monitoring for patients suspected of having life-threatening arrhythmias.

December, 2007



## NorthEast Monitoring, Inc.

Two Clock Tower Place, Suite 555  
Maynard, MA 01754 U.S.A.  
www.nemon.com

phone: 978-461-3992  
fax: 978-461-5991  
toll free: 866-346-5837



0413