

 ■ NEMon.com) (+1) 978-461-3992) Toll Free USA 866-346-5837 	DR300 Holter and Event Recorder with Bluetooth Wireless A revolutionary recorder combining the best of Holter, Event and Mobile Cardiac Telemetry (MCT)	DR200/HE Holter and Event Recorder The industry's smallest, single unit Holter and Event recorder	DR181 3-Channel with OxyHolter Capability Designed for ease-of-use, accuracy of data and durability of the unit	DR181 12-Lead with OxyHolter Capability High-quality research and clinical evaluation tool
⊠ Sales@NEMon.com	CALL TO SERVICE AND ADDRESS OF THE PARTY OF	EVENT		
Current version	4.44	1.11 and 4.44	3.05	3.05
Holter Recorder	X	X	X	X
SD Card data transfer	X	X	X	X
Event Recorder	X	X		
Event, auto-detect	X	X		
3-day 3CH Holter	X	X	X	X
3-day 12-lead			X	X
14-day Holter	X	X		
30-day Event	X	X		
Event with Holter backup	X	X		
Event Bluetooth	X			
Tablet set up interface	X			
2-lead hook-up	Х	X		
3-lead hook-up	X	X		
5-lead hook-up	X	X	X	X
7-lead hook-up	X	X	X	X
12-lead hook-up				X
Oximetry			X	X
Sleep			X	X
Pacemaker detection	X	X	X	X
3-year warranty	X	X	X	X
Batteries	1- AA	1-AA	2-AA	2-AA
Sampling rate	180 samples/second	180 samples/second	1440 samples/second	1440 samples/second

NorthEast Monitoring is a medical electronics company that develops high quality, rugged products to address diagnostic, clinical and research needs in cardiology, respiratory, and pulmonary applications.

Contact NorthEast Monitoring for more information or to place an order.