



**NorthEast Monitoring, Inc.**

*advancing Holter technology™*

## Compare Recorders

**nemon.com**

**(+1) 978-461-3992**

**Toll Free USA 866-346-5837**

**sales@nemon.com**

### **DR400** patch style Holter recorder with Wireless event

Weighing in at just 33 grams, the showerproof 1- or 3-lead patch style recorder solves your remote cardiac monitoring needs



### **DR300** Holter recorder with Wireless event

A revolutionary recorder combining the best of Holter and Event



### **DR200/HE** Holter and Event Recorder

The industry's smallest, single unit Holter and Event recorder



<b>Current version</b>	5.28	4.48	1.11 and 4.47
<b>Holter Recorder</b>	X	X	X
<b>Holter Capacity</b>	14 days	14 days	14 days
<b>Event Recorder</b>	X	X	X
<b>Event, auto-detect</b>	X	X	X
<b>Event capacity</b>	recharge every 14 days	30 days	30 days
<b>Event, wireless</b>	X	X	
<b>Event w/ Holter back up</b>	14 days	14 days	14 days
<b>Pacemaker detection</b>	X	X	X
<b>1-lead hook-up</b>	patch electrode	X	X
<b>3-lead hook-up</b>	patch electrode	X	X
<b>5-lead hook-up</b>		X	X
<b>7-lead hook-up</b>		X	X
<b>12-lead hook-up</b>			
<b>12-lead capacity</b>			
<b>Oximetry</b>			
<b>User interface</b>	PCPatch utility on PC	LCD screen	LCD screen
<b>Data Transfer</b>	Docking station	SD Card	SD Card
<b>Industry-leading warranty</b>	3-year	3-year	3-year
<b>Batteries</b>	Rechargeable internal	1- AA	1-AA
<b>Sampling rate</b>	180 samples/second	180 samples/second	180 samples/second

NorthEast Monitoring is a medical electronics company that develops high quality, rugged products to address diagnostic, clinical and research needs in cardiology, respiratory, and pulmonary applications.

Contact NorthEast Monitoring for more information or to place an order.  
[www.nemon.com](http://www.nemon.com) Updated August 2, 2023 [sales@nemon.com](mailto:sales@nemon.com)